Tierwaage
Animal scale
Pèse-animaux
Báscula de animales
Bilancia di animali
Dierlijke schubben
Łuski zwierząt

6858

Gebrauchsanweisung - Seite 3
Manual - Page 9
Mode d'emploi - Page 15
Instrucciones de uso - Página 21
Istruzioni per l'uso - Pagina 27
Instructies voor gebruik- Pagina 34
Instrukcja użytkowania- Strona 39

DE

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EN

Thank you for purchasing the Soehnle Professional Animal Scale 6858. This scale was developed in accordance with the requirements of practical operation and offers you numerous possibilities to individually adapt the weighing process to your needs in order to make your processes rational and cost-effective.

Please read the operating instructions carefully before use. Improper use could damage you or the terminal.
If you have any questions or have problems with your terminal that are not covered in this manual, please contact your Soehnle Industrial Solutions Service Center or Soehnle Industrial Solutions Customer Service.

Phone: +49 7191 3453-220
Fax: +49 7191 3453-211
E-mail: info@sis.gmbh

Further information and documentation can be found under:

<https://www.soehnle-professional.com/en/site/documents>

Please go to the customer centre on our website www.soehnle-professional.com and select the animal scale 6858 under Downloads / Veterinary scales.

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1. INTRODUCTION

The range of platform scales has a stainless steel top pan on a steel frame and an indicator with a large backlit LCD. The keyboard has 4 easy to use function keys , ,  and . All scales are supplied with the Power Supply module.

2. SPECIFICATIONS

Capacity x Readability d=e=	200 kg x 0.05kg	440 lb x 0.1 lb	7054 oz x 2 oz	440 lb:16 oz x 2 oz
Repeatability	0.05 kg / 0.1 lb			
Linearity	0.10 kg / 0.2 lb			

Platform size	900 x600 x80mm
POWER-supply	4 x1.5V AA battery / DC 5V / USB adaptor 5V (optional)
Indicator Dimensions (w x d x h)	215 x 100x 70 mm
Stabilization time	2-3 seconds
Tare	Full range tare by subtraction
Weighing units	kg, lb, oz, lb.oz
Calibration	Automatic External - user selectable cal weight
Interface	Bi-directional USB
Operating temperature	0°C to 40°C
Humidity	Up to 90% RH non-condensing
Display	25 mm Backlit LCD digits with weight legends for kg, lb, oz, lb.oz and low battery, stable, zero, net weight and Hold symbols
Keypad	Mechanical switches under overlay
Scale housing	Indicator: Aluminium Platform: Mild steel base and stainless steel pan
Applications	Weighing, Dynamic / Animal weighing, Display Hold

3. INSTALLATION

3.1 Unpacking


The scales come with

- | | |
|---|-------------------------------------|
| ✓ | Indicator with bracket and 2 screws |
| ✓ | Weighing platform |
| ✓ | Stainless steel hood |
| ✓ | Rubber mat |
| ✓ | Power supply module |
| ✓ | Instruction manual |

3.2 Locating

- > The scale should not be placed in a location that will reduce the accuracy.
- > Avoid extremes of temperature.
- > Do not place the scale in direct sunlight, open windows, or near air conditioning vents.
- > Avoid unsuitable surfaces. The table or floor should be rigid and free from vibration.
- > Avoid unstable power sources.
- > Do not use the scale near large users of electricity.
- > Do not place the scale near vibrating machinery.
- > Avoid high humidity that might cause condensation.
- > Avoid direct contact with water, e.g. by immersion or spraying.
- > Avoid air movement such as from fans or opening doors.
- > Keep the scales clean.
- > Do not stack material on the scales when they are not in use.

3.3 Setting up


NOTE: When moving the scale, take care to prevent bending of the cable and connector too much as over a period of time this can weaken the wires in the connector and cause malfunction.


4. KEYPAD

	Turns the scale on or off only
	Sends data via USB and combines with Hold functions, if enabled
	Selects the units (kg, lb, lb/oz)
	Sets the display to true zero or net zero by storing the current weight in the tare memory

5. DISPLAY

6. SYMBOLS AND INDICATORS

Symbol	Description
	Scale is in the zero position
	Weighing result is stable
kg, lb, oz or lb.oz	Active weighing unit
	Battery is weak
	Net weight is being displayed
Hold	Display is held as per the Hold parameter setting

7. BATTERY

When the battery needs charging a symbol  on the display will turn on. The battery should be changed when the symbol is on, or connect the power supply module to the rear of the indicator and apply power. The scale does not need to be turned on.

8. BACKLIGHT

The backlight is factory set.

9. AUTO POWER OFF

The auto power off function is factory set.

10. OPERATIONS

10.1 Turning on the scale

- > To turn on the scale, press the  key.
- > The display will show the software revision number and then flash all the digits and symbols before counting down to zero. This ensures all LCD segments are working. The last active weighing unit will be displayed.
- > The scale will turn off automatically to conserve battery.
- > To turn off the scale immediately, press the  key.
- > A battery symbol will be shown when the internal battery needs to be recharged. Connect the power supply module to the rear of the indicator and switch on the mains.

10.2 Zeroing

If the scale is less than 2% of the maximum capacity, pressing  will zero the scale. However, if the weight on the scale is more than 2%, pressing  will tare the scale. See the next section for taring the scale.

10.3 Taring

Press the  key to tare the scale. The weight that was displayed is stored as the tare value. This value is subtracted from the display, leaving zero on the display. The indicator  will be on.

10.4 Selecting the unit

To select the weighing unit press the  key to move from one unit to another. The weighing unit must be enabled by the users beforehand. If a unit is disabled, it cannot be selected by using the  key.

11. ERROR MESSAGES

During the initial power-on testing or during the operation the scale may show an error message. The error messages are described below.

ERROR CODE	DESCRIPTION	POSSIBLE CAUSES
nnnnn	A continuous beep is heard.	Weight on the pan exceeds the capacity of the scale. Remove the weight from the weighing platform.
CALEr	If the selected mass is less than 10% or more than 20% of the capacity of the scale, an error message "CALEr" will be displayed and the scale will return to zero.	Incorrect calibration mass. Repeat the process correctly.

If an error message is shown, repeat the procedure that caused the message such as turning the scale on, calibration or any other functions. If the error message still appears, contact Soehnle Industrial Solutions customer service.

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Soehnle Industrial Solutions GmbH