




Veterinary scale "Easy" with fold function

8320

- ▶ Foldable veterinary scale for small animals
- ▶ Perfect for mobile use
- ▶ Low space requirement, easy to transport
- ▶ Hold and tare function
- ▶ Energy-saving auto-off function



Features

- ▶ Patented folding mechanism with magnetic lock
- ▶ Low space requirement, easy to transport
- ▶ Easy to read LCD display
- ▶ Tare function: Covers/paddings can be tared
- ▶ Hold function: The weight value will be held for a few seconds after unloading
- ▶ Switchable to kg/lb units
- ▶ Energy-saving auto-off function
- ▶ Battery operation with 2 x 1.5 V AAA alkaline batteries (included in scope of delivery)

Products

Order number Variation	Approvable	Division g	Maximum load kg
8320.01.001 Veterinary scale for small animals with fold function	No	10	15

Accessories

Order number
Variation

Description

5040.01.001



Carry bag black, with flexible shoulder strap and two carrying handles
