




Veterinary scale "Easy" with fold function

8320

- ▶ Foldable veterinary scale for small animals
- ▶ Perfect for mobile use
- ▶ Low space requirement, easy to transport
- ▶ Hold and tare function
- ▶ Energy-saving auto-off function



Features

- ▶ Patented folding mechanism with magnetic lock
- ▶ Low space requirement, easy to transport
- ▶ Easy to read LC display
- ▶ Tare function: Covers/paddings can be tared
- ▶ Hold function for easy reading of the weight value after weighing
- ▶ Switchable to kg/lb units
- ▶ Energy-saving auto-off function
- ▶ Battery operation with 2 x 1.5 V AAA alkaline batteries (not included in scope of delivery)

Products

Order number	Approvable	Division	Maximum load
Variation		g	kg
8320.01.001	No	10	15
Veterinary scale for small animals with fold function			

Accessories

Order number

Variation

Description

5040.01.001



Carry bag black, with flexible shoulder strap and two carrying handles
