





Chair weigher

7802

- ▶ Robust construction
- ▶ Simple weighing of persons with a mobility handicap
- ▶ Ergonomically shaped and easy to clean seat
- ▶ Hold and tare function



Features

- ▶ Model 7802: Switchable to kg/lb units
- ▶ Ergonomically shaped and easy to clean seat
- ▶ Dimensions seat: 41 x 40 cm
- ▶ Seat height: approx. 51 cm
- ▶ Hinged arm- and footrests
- ▶ Energy-saving auto-off function
- ▶ Rechargeable battery and mains operation
- ▶ Taring function: Placed underlays can be tared
- ▶ Hold function for easy reading of the weight value after weighing
- ▶ Version .002: two rear swivel castors with rear fixing brakes
- ▶ Version .402: four swivel castors with rear fixing brakes

Products

Order number Variation	Approvable	Division g	Maximum load kg	Model
7802.70.002 Chair weigher with display Basic	No	100	200	2 rear swivel castors with rear fixing brakes
7802.70.402 Chair weigher with display Basic	No	100	200	4 swivel castors with rear fixing brakes