



Digital personal scale

7810

- ▶ Simple to use
- ▶ Tap-on start function
- ▶ Extra flat design
- ▶ Mobile use due to low own weight

Features

The 7810 Digital Personal Scale from Soehnle Professional is characterized by a simple, easy and safe operation. The plastic surface of the weighing platform is very easy to clean. Due to the weighing capacity of up to 150 kg, this digital personal scale is suitable for weighing most of all people. The LCD display with 26 mm high digits enables easy reading of the ascertained weight.

- ▶ Simple to use
- ▶ Instant-start-function by tapping the weighing platform
- ▶ Easy-to-read LC display, digit height 26 mm
- ▶ Units in kg, lb or st, switchable at the bottom of the scale
- ▶ Hold function for easy reading of the weight value after weighing
- ▶ Battery operation with one CR2032 button cell battery (included in scope of delivery)
- ▶ Energy-saving auto-off function

Products

Order number	Approvable	Division	Maximum load	Platform size
Variation		g	kg	mm
7810.01.001	No	100	180	304 x 385 x 26
Digital Personal Scale				
