



Wheelchair weigher

6808

- ▶ Easy to transport
- ▶ Integrated ramps for easy access
- ▶ Functions: Weighing, taring, hold, BMI
- ▶ Transportable due to two castors and carrying handles
- ▶ Energy-saving auto-off function

Features

- ▶ Powder-coated platform with integrated drive-on ramps
- ▶ Passage width: 795 mm
- ▶ Terminal with aluminium housing
- ▶ Large LC display, digit height: 25.4 mm
- ▶ Transportable due to two castors and carrying handles
- ▶ Four leveling feet for levelling and stabilizing the scale
- ▶ Functions: Weighing, taring, hold, BMI
- ▶ Battery level indicator (battery symbol in display)
- ▶ Energy-saving auto-off function
- ▶ Units in kg, lb or lb:oz (factory setting: kg)
- ▶ Serial RS 232 and USB interface (see Accessories for corresponding data cable)
- ▶ Battery operation with four Mignon batteries 1.5 V AA (included in scope of delivery)
- ▶ Power supply operation optional (see Accessories)

Products

Order number	Approvable	Division	Maximum load
Variation		g	kg
6808.01.501	No	100	300
Wheelchir weigher			

Accessories

Order number

Variation

Description

5908.09.001

Power supply 9V / 2A
