



# Veterinary scale

6958

- ▶ Easy to clean surface
- ▶ Weighing platform divided into two parts
- ▶ Functions: Weighing, taring, hold
- ▶ Energy-saving auto-off function
- ▶ Protection class IP65 (for terminal box)

## Features

---

- ▶ Weighing platform divided into two parts
- ▶ Platforms individually packed to reduce shipping costs
- ▶ Four load cells, four adjustable feet and two handles per platform
- ▶ Base frame made of steel
- ▶ Weighing platform made of checker plate
- ▶ Functions: Weighing, taring, hold
- ▶ Battery level indicator (battery symbol in display)
- ▶ Energy-saving auto-off function
- ▶ Units in kg, lb or lb:oz (factory setting: kg)
- ▶ Serial RS 232 and USB interface (see Accessories for corresponding data cable)
- ▶ Battery operation with four Mignon batteries 1.5 V AA (included in scope of delivery)
- ▶ Power supply operation optional (see Accessories)

## Products

---

Order number	Approvable	Division	Maximum load
Variation		g	kg
<b>6958.01.002</b>	No	500	1.500
Veterinary scale			

---

## Accessories

---

**Order number**

Variation

**Description**

---

**5908.09.001**

Power supply 9V / 2A

---

**2550.03.002**

Data cable, 9-pin connection terminal - PC, cable length 2 m

---