



Personal scale with handrail

6803

- ▶ Including support railing for a safe stand
- ▶ Support railing made of stainless steel
- ▶ Functions: Weighing, taring, hold
- ▶ 4 leveling feet for levelling and stabilizing the scale
- ▶ Energy-saving auto-off function

Features

- ▶ Powder-coated base frame, platform made of stainless steel
- ▶ Terminal with aluminium housing
- ▶ Large LC display, digit height: 25.4 mm
- ▶ Stainless steel handrails with padded handles
- ▶ 4 leveling feet for levelling and stabilizing the scale
- ▶ Functions: Weighing, taring, hold
- ▶ Battery level indicator (battery symbol in display)
- ▶ Energy-saving auto-off function
- ▶ Units in kg, lb or lb:oz (factory setting: kg)
- ▶ Serial RS 232 and USB interface (see Accessories for corresponding data cable)
- ▶ Battery operation with four Mignon batteries 1.5 V AA (included in scope of delivery)
- ▶ Power supply operation optional (see Accessories)

Products

Order number Variation	Approvable	Division g	Maximum load kg
6803.01.001 Personal scale with handrail	No	100	300

Accessories

Order number Variation	Description
5908.09.001	Power supply 9V / 600 mA
2557.25.002	UK adapter
2557.25.003	US adapter
2550.03.002	Data cable, 9-pin connection terminal - PC, cable length 2 m