



Personal scale with handrail

6831

- ▶ Including support railing for a safe stand
- ▶ With integrated height rod
- ▶ 4 leveling feet for levelling and stabilizing the scale
- ▶ Functions: Weighing, taring, hold
- ▶ Energy-saving auto-off function

Features

- ▶ Powder-coated platform with integrated support railing
- ▶ Terminal with aluminium housing
- ▶ Large LC display, digit height: 25.4 mm
- ▶ With extendable aluminium height rod, Unit: 1 mm
- ▶ Four leveling feet for levelling and stabilizing the scale
- ▶ Easy transport to the patient due to two castors
- ▶ Functions: Weighing, taring, hold, BMI
- ▶ Battery level indicator (battery symbol in display)
- ▶ Energy-saving auto-off function
- ▶ Units in kg, lb or lb:oz
- ▶ Serial RS 232 and USB interface
- ▶ Battery operation with four Mignon batteries 1.5 V AA (included in scope of delivery)
- ▶ Power supply operation optional (see Accessories)

Products

Order number Variation	Approvable	Division g	Maximum load kg
6831.01.001 Personal scale with handrail	No	100	300

Accessories

Order number	Description
Variation	
5908.09.001	Power supply 9V / 2A
2550.03.014	Data cable RS 232 PC/EDV 3 m