



Veterinary scale

6858

- ▶ For small and medium-sized animals
- ▶ Quick and easy operation
- ▶ Flat stainless steel weighing platform

Features

- ▶ Fast and easy operation via four function keys
- ▶ Easy to read LC display, digit height 25 mm
- ▶ Tare function: Covers/paddings can be tared
- ▶ Hold function for better readability
- ▶ Display unit made of aluminium
- ▶ Large stainless steel weighing platform with separate grey rubber mat 900 x 600 x 80 mm
- ▶ Especially flat
- ▶ Units in kg, lb, oz, lb.oz
- ▶ Energy-saving auto-off function
- ▶ Power supply operation with DC 5V power supply
- ▶ Battery operation optional with four Mignon batteries 1.5 V AA (batteries not included)

Products

Order number	Approvable	Division	Maximum load	Platform size
Variation		g	kg	mm
6858.01.002	No	50	200	900 x 600
Veterinary scale				
