



Baby scale "Cosy"

8310.02

- ▶ Simple operation
- ▶ Easy to read LC display
- ▶ Easy to clean plastic housing
- ▶ Tare, hold and memory function
- ▶ Energy-saving auto-off function



Features

The Soehnle Professional Baby Scale 8310.02 is an elegant baby scale. With its rounded corners and the comfortable plastic bowl, safety and hygiene are top priority. Functions such as the hold function and the memory function for the last weight value enable the baby to be weighed quickly and without stress. In battery mode, the scale can be used on the move and is energy-saving thanks to the integrated auto-off function, as the scale switches off automatically when not in use.

- ▶ Easy to read LC display, digit height 25 mm
- ▶ Tare function: Covers/paddings can be tared
- ▶ Hold function for easy reading of the weight value after weighing
- ▶ Memory function: The last weight value determined can be recalled
- ▶ Easy to clean plastic housing
- ▶ Simple operation
- ▶ Units in kg or lb
- ▶ Energy-saving auto-off function
- ▶ Battery operation with four AA 1.5 V batteries (included)
- ▶ Power supply operation optional (see accessories)

Products

Order number	Approvable	Division	Maximum load
Variation		g	kg
8310.02.001 Baby scale "Cosy"	No	5	20

Accessories

Order number

Variation

Description

5908.09.002

Power supply 6V / 500 mA
