



Chair weigher conformity assessed

6702

- ▶ Conformity assessed according to class III, MPG class 1 m
- ▶ Ergonomically shaped and easy to clean seat
- ▶ Hinged arm- and footrests
- ▶ BMI function
- ▶ Four heavy-duty swivel castors with parking brakes front and rear



Features

- ▶ Conformity assessed according to class III, MPG class 1 m
- ▶ Ergonomically shaped and easy to clean seat
- ▶ Dimensions seat: 430 x 420 x 80 mm
- ▶ Seat height: approx. 54 cm
- ▶ Hinged arm- and footrests
- ▶ Large LC display with backlight, digit height: 20 mm
- ▶ Energy-saving auto-off function
- ▶ Battery and mains operation 100 - 240 V, 50-60 Hz (included in scope of delivery)
- ▶ Taring function: Placed underlays can be tared
- ▶ Hold function for easy reading of the weight value after weighing
- ▶ BMI function: Automatic calculation of the BMI by entering the body height and the determined weight value
- ▶ Four heavy-duty swivel castors with parking brakes front and rear
- ▶ Units in kg

Products

Order number	Approvable	Division	Maximum load
Variation		g	kg
6702.01.402	Yes	100	300
Chair weigher			
