



Wheelchair weigher Conformity assigned

6708

- ▶ Conformity assessed according to Class III, MPG Class 1 m
- ▶ Platform with integrated drive-on ramps
- ▶ Including foldable railing, white painted
- ▶ Transportable due to two integrated castors and handle
- ▶ Battery and mains operation



Features

- ▶ Conformity assessed according to Class III, MPG Class 1 m
- ▶ Platform with integrated drive-on ramps
- ▶ Including foldable railing, white painted
- ▶ Passage width: 760 mm
- ▶ Terminal with plastic housing
- ▶ Backlit LC display, digit height: 20 mm
- ▶ Transportable due to two integrated castors and handle
- ▶ Functions: weighing, taring, hold, BMI
- ▶ Battery level indicator (battery symbol in the display)
- ▶ Energy-saving auto-off function
- ▶ Units in kg
- ▶ Battery and mains operation 100 - 240 V, 50-60 Hz (included in scope of delivery)

Products

Order number	Approvable	Division	Maximum load	Platform size	Product weight
Variation		g	kg	mm	kg
6708 Wheelchair weigher	Yes	100	300	970x840x82	46