



# Mechanical personal scale

6163

- ▶ Classic circle pointer display
- ▶ Easy to read in kg-steps

## Features

---

- ▶ Large and clear circle pointer display, diameter: 18 cm
- ▶ Easy to read in kg-steps
- ▶ Manual zero resetting possible
- ▶ Attractive chrome trim ring, diameter: 21 cm
- ▶ Sturdy steel housing
- ▶ Weighing platform with easy-to-clean plastic covering

## Products

---

Order number	Approvable	Division	Maximum load
Variation		g	kg
<b>6163.02.001</b>	No	1000	160
Mechanical personal scale			

---